


Description

Product Name	YL-109
Brief Description	Inhibitors
Purification	99.77%
Target Name	AhR activator
Calculated MW	257.31
Formulation	C ₁₄ H ₁₁ NO ₂ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: COC1=C/C(=C/2\NC3=CC=CC=C3S2)/C=CC1=O
 Pathways: Metabolism
 Receptor: AhR
 Boiling pt: Melting
 pt: Solubility: Soluble in Water
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Hiyoshi H, et al. Sci Rep. 2014 Nov 18;4:7095.

Note: This product is for in vitro research use only