

## INCB28060(Capmatinib)

Catalog No: #S1963


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	INCB28060(Capmatinib)
Brief Description	Inhibitors
Purification	98.00%
Target Name	c-Met/HGFR inhibitor
Calculated MW	412.42
Formulation	C23H17FN6O
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Other  
 SMILES:CNC(=O)C1=C(F)C=C(C=C1)C1=NN2C(CC3=CC4=C(C=C3)N=CC=C4)=CN=C2N=C1  
 Pathways:Tyrosine Kinase/AdaptorsReceptor:c-Met  
 Boling pt:Melting pt:Solubility:DMSO: 2 mg/mL (4.84 mM)  
 Appearance:Yellow Powder  
 Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.  
 Stock solution can be stored below -20 °C for several month.

## References

1. Liu X, et al, Clin Y Res, 2011, 17(227), 7127-7138.

Note: This product is for in vitro research use only