

Bikinin

Catalog No: #S1965


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bikinin
Brief Description	Inhibitors
Purification	99.50%
Target Name	GSK-3 inhibitor; Others inhibitor
Calculated MW	273.08
Formulation	C9H9BrN2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: OC(=O)CCC(=O)NC1=NC=C(Br)C=C1
 Pathways: PI3K/Akt/mTOR signaling; Others; Stem Cells
 Receptor: GSK-3; Shaggy-like kinase
 Boiling pt: Melting pt: Solubility: DMSO: 54 mg/mL (197.74 mM)
 Appearance: White Cotton Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. De Rybel B, et al. Chem Biol. 2009, 16(6), 594-604.

Note: This product is for in vitro research use only