

DBEQ


Catalog No: #S1969

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DBEQ
Brief Description	Inhibitors
Purification	99.15%
Target Name	p97 inhibitor
Calculated MW	340.42
Formulation	C22H20N4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C1=NC2=CC=CC=C2C(NCC2=CC=CC=C2)=N1C1=CC=CC=C1
Pathways: Ubiquitination Receptor: p97
Boiling pt: Melting pt: Solubility: Ethanol: 5 mg/mL (14.68 mM); DMSO: 68 mg/mL (199.75 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Chou TF, et al. Proc Natl Acad Sci U S A, 2011, 108(12), 4834-4839.

Note: This product is for in vitro research use only