


Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WAY-262611
Brief Description	Inhibitors
Purification	98.00%
Target Name	Wnt/beta-catenin agonist
Calculated MW	318.42
Formulation	C20H22N4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: NCC1CCN(CC1)C1=NC=CC(=N1)C1=CC2=CC=CC=C2C1
Pathways: Cytoskeletal Signaling; Stem Cells
Receptor: Wnt/ β -catenin
Boiling pt: Melting pt: Solubility: Soluble in DMSO
Appearance: Off-white Color Solid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Pelletier JC, et al. J Med Chem. 2009 Nov 26;52(22):6962-5.

Note: This product is for in vitro research use only