


## Description

Product Name	PFI-4
Brief Description	Inhibitors
Purification	98.00%
Target Name	Epigenetic Reader Domain inhibitor
Calculated MW	380.44
Formulation	C21H24N4O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: COC1=CC=CC=C1C(=O)NC1=C(C=C2N(C)C(=O)N(C)C2=C1)N1CCCC1 Pathways: Chromatin/Epigenetic Receptor: BRPF1; BRPF2; BRPF3  
 Boiling pt: Melting pt: Solubility: DMSO: 50 mM Appearance: Faint Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

- Demont EH, et al. ACS Med Chem Lett. 5(11):1190-1195.

Note: This product is for in vitro research use only