


## Description

Product Name	GW 627368X
Brief Description	Inhibitors
Purification	99.92%
Target Name	Prostaglandin Receptor antagonist
Calculated MW	544.62
Formulation	C30H28N2O6S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: CCOC1=C2CN(C(=O)C2=C(OCC)C2=CC=CC=C12)C1=CC=C(CC(=O)NS(=O)(=O)C2=CC=CC=C2)C=C1 Pathways: GPCR/G Protein Receptor: EP4 Boiling pt: Melting pt: Solubility: DMSO: 10 mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Wilson RJ, et al. Br J Pharmacol. 2006 Jun; 148(3): 326-39.

Note: This product is for in vitro research use only