

Cinaciguat

Catalog No: #S1984


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cinaciguat
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	565.7
Formulation	C36H39NO5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: OC(=O)CCCCN(CCC1=CC=CC=C1OCC1=CC=C(CCC2=CC=CC=C2)C=C1)CC1=CC=C(C=C1)C(O)=O Pathways: Others Recept
or: sGC Boling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Mihai Gheorghiade¹, et al. Eur J Heart Fail (2012) 14 (9): 1056-1066.

Note: This product is for in vitro research use only