

## Luminespib (AUY-922, NVP-AUY922)

Catalog No: #S1989


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Luminespib (AUY-922, NVP-AUY922)
Brief Description	Inhibitors
Purification	98.97%
Target Name	HSP inhibitor
Calculated MW	465.54
Formulation	C26H31N3O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CCNC(=O)C1=C(/C(=C/C2=CC(=O)O)C(C)C)/ON1)C3=CC=C(C=C3)CN4CCOCC4  
 Pathways: Cytoskeletal Signaling; Metabolism  
 Receptor: HSP90α; HSP90β  
 Boiling pt: Melting pt: Solubility: Ethanol: 31 mg/mL (66.58 mM); DMSO: 93 mg/mL warmed (199.76 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Wainberg ZA, et al. Mol Y Ther. 2013 Apr;12(4):509-19.

Note: This product is for in vitro research use only