


Description

Product Name	JNJ-31020028
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others antagonist
Calculated MW	565.68
Formulation	C34H36FN5O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CCN(CC)C(=O)C(N1CCN(CC1)C1=C(F)C=C(NC(=O)C2=CC=CC=C2C2=CN=CC=C2)C=C1)C1=CC=CC=C1 Pathways: Others Receptor: Y2 receptor Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Shoblock JR, et al. *Psychopharmacology (Berl)*. 2010 Feb;208(2):265-77.

Note: This product is for in vitro research use only