

## Ciproxifan maleate

Catalog No: #S2004


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ciproxifan maleate
Brief Description	Inhibitors
Purification	99.45%
Target Name	Histamine Receptor antagonist
Calculated MW	386.4
Formulation	C <sub>20</sub> H <sub>22</sub> N <sub>2</sub> O <sub>6</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: C1CC1C(=O)C2=CC=C(C=C2)OCCCC3=CN=CN3.C(=C\C(=O)O)\C(=O)O  
 Pathways: GPCR/G Protein;  
 Immunology/Inflammation; Neuroscience  
 Receptor: H3 receptor  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Motawaj M, Arrang JM. Eur J Neurosci. 2011 Apr;33(7):1197-204.

Note: This product is for in vitro research use only