

Fexaramine


Catalog No: #S2015

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fexaramine
Brief Description	Inhibitors
Purification	98.00%
Target Name	FXR agonist
Calculated MW	496.64
Formulation	C32H36N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: OtherSMILES: CN(C)C1=CC=C(C=C1)C2=CC=C(C=C2)CN(C3=CC=CC(=C3)/C=C/C(=O)OC)C(=O)C4CCCCC4Pathways: MetabolismReceptor: FXRBoling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Lam IP, et al. Am J Physiol Gastrointest Liver Physiol. 2009 Jul;297(1): G90-7.

Note: This product is for in vitro research use only