


Description

Product Name	Apilimod
Brief Description	Inhibitors
Purification	98.00%
Target Name	IL Receptor inhibitor
Calculated MW	418.49
Formulation	C23H26N6O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC1=CC=CC(=C1)/C=N/NC2=NC(=NC(=C2)N3CCOCC3)OCC4=CC=CC=N4 Pathways: Immunology/Inflammation Receptor: IL-12; IL-23 Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Off-white Color Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Wada Y, et al. PLoS One. 2012;7(4):e35069.

Note: This product is for in vitro research use only