

Dihydroethidium

Catalog No: #S2021


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dihydroethidium
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	315.41
Formulation	C ₂₁ H ₂₁ N ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CCN1C(C2=CC=CC=C2)C2=C(C=CC(N)=C2)C2=C1C=C(N)C=C2
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Purple Solid
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zielonka J, et al. J Biol Chem. 2012 Jan 27;287(5):2984-95.

Note: This product is for in vitro research use only