


Description

Product Name	CZC-54252
Brief Description	Inhibitors
Purification	98.00%
Target Name	LRRK2 inhibitor
Calculated MW	504.39
Formulation	C22H25ClN6O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other
 SMILES: Cl.COC1=C(NC2=NC=C(C1)C(NC3=CC=CC=C3NS(C)(=O)=O)=N2)C=CC(=C1)N1CCOCC1
 Pathways: Autophagy Receptor: LRRK2
 Boiling pt: Melting pt: Solubility: DMSO: 100 mM
 Appearance: Faint Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Liu F, et al. J Pharmacol Exp Ther, 2008, 327(3), 827-839.

Note: This product is for in vitro research use only