


Description

Product Name	CYM5442
Brief Description	Inhibitors
Purification	98.00%
Target Name	S1P Receptor agonist
Calculated MW	409.49
Formulation	C23H27N3O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:Cl.CCOC1=C(OCC)C=C(C=C1)C1=NC(=NO1)C1=C2CCC(NCCO)C2=CC=C1Pathways:GPCR/G

ProteinReceptor:SphKBoling pt:Melting pt:Solubility:DMSO: > 10 mMAppearance:Grayish Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Poti F, et al. Vascul Pharmacol. 2012 Aug 19;57(1):56-64.

Note: This product is for in vitro research use only