

Bohemine

Catalog No: #S2029


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Bohemine |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | CDK inhibitor |
| Calculated MW | 340.42 |
| Formulation | C18H24N6O |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Other
 SMILES: CC(C)N1C=NC2=C1N=C(NCCCO)N=C2NCC1=CC=CC=C1
 Pathways: Cell Cycle/Checkpoint Receptor: CDK
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Liskova B, et al. Chem Res Toxicol. 2012 Feb 20;25(2):500-9.

Note: This product is for in vitro research use only