

Tiplaxtinin(PAI-039)

Catalog No: #S2030


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tiplaxtinin(PAI-039)
Brief Description	Inhibitors
Purification	99.00%
Target Name	PAI-1 inhibitor
Calculated MW	439.38
Formulation	C ₂₄ H ₁₆ F ₃ NO ₄
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: OC(=O)C(=O)C1=CN(CC2=CC=CC=C2)C2=C1C=C(C=C2)C1=CC=C(OC(F)(F)F)C=C1 Pathways: Metabolism Receptor: PAI-1 Boli
 ng pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Elokda H, et al. J Med Chem. 2004 Jul 1;47(14):3491-4.

Note: This product is for in vitro research use only