


Description

Product Name	AG-1557
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	407.21
Formulation	C16H14IN3O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:COC1=C(OC)C=C2C(NC3=CC(I)=CC=C3)=NC=NC2=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Faint Yellow SolidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only