


Description

Product Name	GTPL5846
Brief Description	Inhibitors
Purification	98.00%
Target Name	GPR agonist
Calculated MW	239.31
Formulation	C12H21N3O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CCCCCCCNC1=CC(=O)NC(=O)N1
Pathways: Endocrinology/Hormones
Receptor: GPR84
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Faint Yellow Solid
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Suzuki M, et al. J Biol Chem. 2013 Apr 12;288(15): 10684-91.

Note: This product is for in vitro research use only