


Description

Product Name	IDH1-IN-1
Brief Description	Inhibitors
Purification	99.50%
Target Name	Dehydrogenase inhibitor
Calculated MW	498.59
Formulation	C30H31FN4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CC1=CC=CC=C1C(N(C(=O)CN1C=NC2=CC=CC=C2)C1=CC(F)=CC=C1)C(=O)NC1CCCC1 Pathways: Metabolism Receptor: mutant IDH1 R132H
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Crystalline Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Janeta Popovici-Muller, et al. Patent WO2012009678 A1.

Note: This product is for in vitro research use only