

Varespladib (LY315920)

Catalog No: #S2044


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Varespladib (LY315920)
Brief Description	Inhibitors
Purification	99.70%
Target Name	Phospholipase inhibitor
Calculated MW	380.39
Formulation	C21H20N2O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

System SMILES: CCC1=C(C(=O)C(N)=O)C2=C(C=CC=C2OCC(O)=O)N1CC1=CC=CC=C1 Pathways: Metabolism Receptor: hnsPLA2 Boiling pt: Melting pt: Solubility: DMSO: 76 mg/mL (199.79 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Draheim SE, et al. J Med Chem. 1996, 39(26), 5159-5175.

Note: This product is for in vitro research use only