


## Description

Product Name	UNC0064-12
Brief Description	Inhibitors
Purification	97.00%
Target Name	VEGFR inhibitor
Calculated MW	364.45
Formulation	C19H24N8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer  
SMILES: NCCCNC1=CC=C(NC2=NC=CC(NC3=NNC(=C3)C3CC3)=N2)C=C1  
Pathways: Angiogenesis; Tyrosine Kinase/Adaptors Receptor: VEGFR  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. MULTIPLEXED KINASE INHIBITOR BEADS AND USES THEREOF. Patent Application WO/2013/055780 Kind Code: A1

Note: This product is for in vitro research use only