


Description

Product Name	ACT 058362; Palosuran
Brief Description	Inhibitors
Purification	98.04%
Target Name	GPR antagonist
Calculated MW	418.53
Formulation	C25H30N4O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CC1=NC2=CC=CC=C2C(NC(=O)NCCN2CCC(O)(CC3=CC=CC=C3)CC2)=C1 Pathways: Endocrinology/Hormones Receptor: Urotensin Receptor (GPR14) Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Clozel M, et al. J Pharmacol Exp Ther. 2004 Oct;311(1):204-12.

Note: This product is for in vitro research use only