

## Quizartinib (AC220)

Catalog No: #S2066


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Quizartinib (AC220)
Brief Description	Inhibitors
Purification	99.42%
Target Name	FLT inhibitor
Calculated MW	560.67
Formulation	C29H32N6O4S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Cancer  
 SMILES: CC(C)(C)C1=CC(NC(=O)NC2=CC=C(C=C2)C2=CN3C(SC4=C3C=CC(OCCN3CCOCC3)=C4)=N2)=NO1  
 Pathways: Angiogenesis; Tyrosine Kinase/Adaptors/Receptor: FLT3 (ITD); FLT3 (WT)  
 Boiling pt: Melting pt: Solubility: DMSO: 33.2 mg/mL (59.21 mM)  
 Appearance: Yellow Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Zarrinkar PP, et al. Blood, 2009, 114(14), 2984-2992.

Note: This product is for in vitro research use only