

MCOPPB trihydrochloride

Catalog No: #S2071


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	MCOPPB trihydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Opioid Receptor agonist
Calculated MW	518.01
Formulation	C ₂₆ H ₄₃ Cl ₃ N ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CC1(CCCCCC1)N2CCC(CC2)N3C4=CC=CC=C4N=C3[C@@H]5CCCCNC5.Cl.Cl.Cl Pathways: Endocrinology/Hormones; GPCR/G Protein; Neuroscience Receptor: Opioid Receptor Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Hayashi, et al. Journal of Medicinal Chemistry (2009), 52(3), 610-625.

Note: This product is for in vitro research use only