

## B-Raf inhibitor 1

Catalog No: #S2074


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	B-Raf inhibitor 1
Brief Description	Inhibitors
Purification	98.78%
Target Name	Raf inhibitor
Calculated MW	478.94
Formulation	C26H19ClN8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: CC1=C(NC2=C(C=CC=N2)C2=C3NC=NC3=NC=N2)C2=C(C=C1)C(NC1=CC=C(Cl)C=C1)=NC=C2 Pathways: MAPK Receptor: B-Raf  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Wang X, et al. J Med Chem. 2012 Sep 13;55(17): 7332-41.

Note: This product is for in vitro research use only