

Pimavanserin

Catalog No: #S2076


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pimavanserin
Brief Description	Inhibitors
Purification	98.00%
Target Name	5-HT Receptor agonist
Calculated MW	427.55
Formulation	C25H34FN3O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous systemSMILES:CC(C)COC1=CC=C(CNC(=O)N(CC2=CC=C(F)C=C2)C2CCN(C)CC2)C=C1Pathways:GPCR/G Protein;
NeuroscienceReceptor:5-HT2ABoling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility ,
please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Vanover KE, et al. J Pharmacol Exp Ther. 2006 May;317(2):910-8.

Note: This product is for in vitro research use only