

## Eniporide

Catalog No: #S2089


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Eniporide
Brief Description	Inhibitors
Purification	98.00%
Target Name	ATPase inhibitor
Calculated MW	320.37
Formulation	C <sub>14</sub> H <sub>16</sub> N <sub>4</sub> O <sub>3</sub> S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CC1=CC(N2C=CC=C2)=C(C=C1C(=O)N=C(N)N)S(=O)(=O)C  
 Pathways: Membrane transporter/Ion channel  
 Receptor: Na<sup>+</sup>/K<sup>+</sup>-ATPase  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
 Appearance: Pink Solid  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Zeymer U, et al. J Am Coll Cardiol. 2001 Nov 15;38(6):1644-50.

Note: This product is for in vitro research use only