

Andarine

Catalog No: #S2108


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Andarine
Brief Description	Inhibitors
Purification	99.84%
Target Name	Androgen Receptor agonist
Calculated MW	441.36
Formulation	C19H18F3N3O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Nervous

systemSMILES:CC(=O)NC1=CC=C(C=C1)OC[C@@](C)(C(=O)NC2=CC(=C(C=C2)[N+](=O)[O-])C(F)(F)F)O

Pathways:Endocrinology/HormonesRece

ptor:Androgen ReceptorBoiling pt:Melting pt:Solubility:Ethanol: 88 mg/mL (199.38 mM); DMSO: 88 mg/mL (199.38 mM)Appearance:Yellow

PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be

stored below -20 °C for several month.

References

1. Yin D, et al. J Pharmacol Exp Ther. 2003, 304(3), 1334-1340.

Note: This product is for in vitro research use only