

Fexinidazole

Catalog No: #S2112


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fexinidazole
Brief Description	Inhibitors
Purification	99.79%
Target Name	Others
Calculated MW	279.31
Formulation	C ₁₂ H ₁₃ N ₃ O ₃ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: CSC1=CC=C(OCC2=NC=C(N2C)[N+](=O)[O-])C=C1
Pathways: Others
Receptor: antiparasitic
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Raether W, Seidenath H. Ann Trop Med Parasitol. 1983 Feb;77(1):13-26.

Note: This product is for in vitro research use only