

Trametinib (GSK1120212)

Catalog No: #S2125


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Trametinib (GSK1120212)
Brief Description	Inhibitors
Purification	99.39%
Target Name	MEK inhibitor
Calculated MW	615.39
Formulation	C26H23FIN5O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: CN1C(=O)C(C)=C2N(C(=O)N(C3CC3)C(=O)C2=C1NC1=C(F)C=C(I)C=C1)C1=CC(NC(C)=O)=CC=C1 Pathways: MAPK Receptor: MEK1; MEK2 Boiling pt: Melting pt: Solubility: DMSO: 22 mg/mL warmed (35.74 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Yamaguchi T, et al. Int J Oncol, 2011, 39(1), 23-31.

Note: This product is for in vitro research use only