


## Description

Product Name	SBI-0206965
Brief Description	Inhibitors
Purification	98.00%
Target Name	Autophagy antagonist
Calculated MW	489.32
Formulation	C21H21BrN4O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: OtherSMILES:CCCCCCCC=CCCCCCCC(C(=O)OCC(C)(COC(=O)C(CCCCCC=CCCCCCCC)NC(=O)OC(C)(C)C)OC(=O)C(CCCCCC=CCCCCCCC)NC(=O)OC(C)(C)C)NC(=O)OC(C)(C)CPathways:AutophagyReceptor:ULK1Boiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:Off-white Color Spongiform SolidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Egan DF, et al. Mol Cell. 2015 Jul 16;59(2):285-97.

Note: This product is for in vitro research use only