

Zaprinast

Catalog No: #S2129


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Zaprinast
Brief Description	Inhibitors
Purification	98.00%
Target Name	PDE inhibitor
Calculated MW	271.27
Formulation	C13H13N5O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory systemSMILES:CCOC1=CC=CC=C1C1=NC(=O)C2=NNNC2=N1Pathways:MetabolismReceptor:PDEBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Taniguchi Y, et al. FEBS Lett. 2006 Sep 18;580(21):5003-8.

Note: This product is for in vitro research use only