

Propidium Iodide

Catalog No: #S2130


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Propidium Iodide
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	668.39
Formulation	C27H34I2N4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:OtherSMILES:[I-].[I-].CC[N+](C)(CC)CCC[N+]=C2=C(N)C=CC2=CC(N)=CC2=C1C1=CC=CC=C1Pathways:OthersReceptor:OthersBolin
g pt:Melting pt:Solubility:Water: 5 mM (sonication)Appearance:Red PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C
and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Uchikura K, et al. Zygote. 2012 Nov;20(4):333-7.

Note: This product is for in vitro research use only