

Geniposidic Acid

Catalog No: #S2136


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Geniposidic Acid
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	374.34
Formulation	C16H22O10
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: C1C=C([C@@H]2[C@H]1C(=CO[C@H]2O[C@H]3[C@@H]([C@H]([C@@H]([C@H](O3)CO)O)O)C(=O)O)COP

Pathways: Others Receptor: Others Boling pt: Melting pt: Solubility: Ethanol: 2 mg/mL (5.34 mM); DMSO: 16 mg/mL (42.74 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Kim SJ, et al. J Ethnopharmacol. 2013 Mar 7;146(1):271-7.

Note: This product is for in vitro research use only