

Butacaine

Catalog No: #S2139


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Butacaine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Sodium Channel inhibitor
Calculated MW	306.44
Formulation	C18H30N2O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CCCCN(CCCC)CCCOC(=O)C1=CC=C(C=C1)NPathways:Membrane Transporter/Ion ChannelReceptor:Sodium ChannelBoiling pt:Melting pt:Solubility:Soluble in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Avdeef A, et al. Eur J Pharm Sci. 2008 Jan;33(1):29-41.

Note: This product is for in vitro research use only