

Mecamylamine hydrochloride


Catalog No: #S2141

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mecamylamine hydrochloride
Brief Description	Inhibitors
Purification	99.50%
Target Name	AChR antagonist
Calculated MW	203.75
Formulation	C11H22ClN
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: Cl.CNC1(C)C2CCC(C2)C1(C)C
Pathways: Neuroscience
Receptor: nAChR
Boiling pt: Melting pt:
Solubility: Soluble in DMSO
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Struthers AM, et al. *Pharmacol Biochem Behav.* 2009 Dec;94(2):319-28.

Note: This product is for in vitro research use only