

## Bergapten

Catalog No: #S2143


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Bergapten
Brief Description	Inhibitors
Purification	99.95%
Target Name	P450 inhibitor
Calculated MW	216.19
Formulation	C12H8O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer SMILES: COC1=C2C=CC(=O)OC2=CC=C1C=CO Pathways: Metabolism Receptor: CYP Boling pt: Melting pt: Solubility: DMSO: 5 mg/mL (23.12 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Panno ML, et al. Breast Y Res Treat, 2012, 136(2), 443-455.

Note: This product is for in vitro research use only