

Fulvestrant

Catalog No: #S2146


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fulvestrant
Brief Description	Inhibitors
Purification	99.89%
Target Name	Estrogen/progestogen Receptor antagonist
Calculated MW	606.77
Formulation	C32H47F5O3S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@]12CC[C@H]3[C@H]([C@@H]1CC[C@@H]2O)[C@@H](CC4=C3C=CC(=C4)O)CCCCCCCCCS(=O)CCCC(C(F)(F)F)(F)

FPathways:Endocrinology/HormonesReceptor:ERBoiling pt:Melting pt:Solubility:Ethanol: 100 mg/mL (164.8 mM); DMSO: 100 mg/mL (164.8

mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wakeling AE, et al. Y Res, 1991, 51(15), 3867-3873.

Note: This product is for in vitro research use only