

Nutlin-3


Catalog No: #S2158

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nutlin-3
Brief Description	Inhibitors
Purification	99.89%
Target Name	Mdm2 inhibitor
Calculated MW	581.49
Formulation	C30H30Cl2N4O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other
 SMILES: COC1=CC(OC(C)C)=C(C=C1)C1=NC(C(N1C(=O)N1CCNC(=O)C1)C1=CC=C(Cl)C=C1)C1=CC=C(Cl)C=C1
 Pathways: Apoptosis
 Receptor: p53-MDM2
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Stéphane Supiot, et al. *Molecular Therapy* (2008), 7(4), 993-999.

Note: This product is for in vitro research use only