

CaMKP Inhibitor

Catalog No: #S2159


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CaMKP Inhibitor
Brief Description	Inhibitors
Purification	99.98%
Target Name	CaMK inhibitor
Calculated MW	377.32
Formulation	C10H12NNaO9S2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: [Na+].NC1=C(C=C(C2=C1C(O)=CC=C2)S([O-])(=O)=O)S(O)(=O)=O
 Pathways: Neuroscience
 Receptor: CaMKP
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Gray Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only