

## Cardiogenol C hydrochloride

Catalog No: #S2161


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Cardiogenol C hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	296.75
Formulation	C13H17ClN4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: Cl.COC1=CC=C(NC2=NC=CC(NCCO)=N2)C=C1  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting  
 pt: Solubility: DMSO: 10 mM  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Mike AK, et al. Cell Physiol Biochem. 2014;33(1):205-21.

Note: This product is for in vitro research use only