

Maribavir

Catalog No: #S2162


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Maribavir
Brief Description	Inhibitors
Purification	99.85%
Target Name	HCV Protease inhibitor
Calculated MW	376.24
Formulation	C ₁₅ H ₁₉ Cl ₂ N ₃ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC(C)NC1=NC2=CC(=C(C=C2N1[C@@H]3[C@H]([C@H]([C@@H](O3)CO)O)O)Cl)Cl Pathways: Microbiology & virology; Proteases/Proteasome Receptor: HCMV Boiling pt: 611°C Melting pt: Solubility: DMSO: 50mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Karen KB, et al. Antimicrob Agents Chem Other. August 2002 vol. 46 no. 8 2365-2372
2. Whitehurst CB, et al. MJ Virol. 2013 May; 87(9):5311-5.

Note: This product is for in vitro research use only