

Dihydrocapsaicin

Catalog No: #S2163


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dihydrocapsaicin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	307.434
Formulation	C18H29NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OthersSMILES:CC(C)CCCCCCC(=O)NCC1=CC(=C(C=C1)O)OCPathways:OthersReceptor:VR1Boiling pt:Melting pt:Solubility:Water: limited solubleAppearance:Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lawson T, et al. The mutagenicity of capsaicin and dihydrocapsaicin in V79 cells. Cancer Lett. 1989 Nov 30;48(2):109-13.

Note: This product is for in vitro research use only