

Wiskostatin

Catalog No: #S2164


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Wiskostatin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others inhibitor
Calculated MW	426.15
Formulation	C ₁₇ H ₁₈ Br ₂ N ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: CN(C)CC(O)CN1C2=C(C=C(Br)C=C2)C2=C1C=CC(Br)=C2
 Pathways: Others
 Receptor: WASP(Arp)2/3
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Guerriero CJ, Weisz OA. Am J Physiol Cell Physiol. 2007 Apr;292(4):C1562-6.

Note: This product is for in vitro research use only