

## SKF-96365 hydrochloride

Catalog No: #S2170


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	SKF-96365 hydrochloride
Brief Description	Inhibitors
Purification	99.48%
Target Name	TRP/TRPV Channel inhibitor
Calculated MW	402.91
Formulation	C22H27ClN2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: COC1=CC=C(C=C1)CCCOC(CN2C=CN=C2)C3=CC=C(C=C3)OC. CI Pathways: Membrane transporter/Ion channel Receptor: TRP Boiling pt: Melting pt: Solubility: Water: 100 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Singh A, et al. Br J Pharmacol. 2010, 160(6):1464-1475.

Note: This product is for in vitro research use only