


## Description

Product Name	SQ 22536
Brief Description	Inhibitors
Purification	98.80%
Target Name	AChR inhibitor
Calculated MW	205.22
Formulation	C9H11N5O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Others SMILES: C1CC(OC1)N2C=NC3=C2N=CN=C3N Pathways: Neuroscience Receptor: adenylyate cyclase Boiling pt: Melting pt: Solubility: DMSO > 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Fabbri E, et al. J Enzyme Inhib. 1991;5(2):87-98.

Note: This product is for in vitro research use only