

## Apigenin

Catalog No: #S2175


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Apigenin
Brief Description	Inhibitors
Purification	97.95%
Target Name	P450 inhibitor
Calculated MW	270.24
Formulation	C15H10O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Metabolism systemSMILES:OC1=CC=C(C=C1)C1=CC(=O)C2=C(O)C=C(O)C=C2O1Pathways:MetabolismReceptor:CYP2C9Boiling pt:Melting pt:Solubility:DMSO: 54 mg/mL (199.82 mM)Appearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Lin JK, et al. J Cell Biochem Suppl, 1997, 28-29, 39-48.

Note: This product is for in vitro research use only