

Kaempferol

Catalog No: #S2177


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Kaempferol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Estrogen/progestogen Receptor antagonist
Calculated MW	286.23
Formulation	C15H10O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory

systemSMILES:OC1=CC=C(C=C1)C1=C(O)C(=O)C2=C(O)C=C(O)C=C2O1Pathways:Endocrinology/HormonesReceptor:ERBoiling pt:Melting pt:Solubility:Ethanol: 57 mg/mL (199.14 mM); DMSO: 57 mg/mL (199.14 mM)Appearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lin F, et al. Int Immunopharmacol. 2015 Oct;28(2):859-65.

Note: This product is for in vitro research use only